

MAVAND Solutions GmbH
 Ulrichstr. 23
 72116 Mössingen
 Tel.: 07473 / 9580 - 28
 Fax: 07473 / 9580 - 29

Measurement
 Rapid STAT
 New Collector with
 5 different cannabis user

2007 

User 1 dated 24.02.2007				
Nr.	Time for collecting saliva after smoking	New Rapid STAT Collector Mirco flocced swab		
		collected quantity of saliva ml	result Rapid STAT THC	Result GC/MS THC ng/ml
4	2 hour	0,22	pos	49,3
5	3 hour	0,18	pos	39,9
6	4 hour	0,18	pos	39,4
7	5 hour	0,20	pos	27,8
8	6 hour	0,24	pos	25,1
9	7 hour	0,21	pos	17,3

User 2 dated 10.03.2007				
Nr.	Time for collecting saliva after smoking	New Rapid STAT Collector Mirco flocced swab		
		collected quantity of saliva ml	result Rapid STAT THC	Result GC/MS THC ng/ml
2	30 min.	0,21	pos	375,0
3	1 hour	0,19	pos	18,6
4	2 hour	0,14	pos	13,5
5	3 hour	0,26	pos	8,6
6	4 hour	0,23	pos	15,2
7	5 hour	0,25	pos	9,4
8	6 hour	0,24	pos	8,9
9	7 hour	0,21	pos	7,9
10	8 hour	0,24	pos	7,4

User 3 dated 16.06.2007				
Nr.	Time for collecting saliva after smoking	New Rapid STAT Collector Mirco flocced swab		
		quantity of saliva ml	result Rapid STAT THC	Result GC/MS THC ng/ml
2	30 min.	0,30	pos	445,0
3	1 hour	0,25	pos	285,0
4	2 hour	0,19	pos	67,4
5	3 hour	0,21	pos	33,9
6	4 hour	0,24	pos	10,8
7	5 hour	0,27	pos	11,4
8	6 hour	0,28	pos	12,4
9	7 hour	0,30	pos	14,9
10	8 hour	0,28	pos	5,3
11	9 hour	0,27	pos	4,5

MAVAND Solutions GmbH
 Ulrichstr. 23
 72116 Mössingen
 Tel.: 07473 / 9580 - 28
 Fax: 07473 / 9580 - 29

Measurement
 Rapid STAT
 New Collector with
 5 different cannabis user

2007



User 4 dated 17.06.2007				
Nr.	Time for collecting saliva after smoking	New Rapid STAT Collector Mirco flocced swab		
		quantity of saliva ml	result Rapid STAT THC	Result GC/MS THC ng/ml
2	30 min.	0,28	pos	298,3
3	1 hour	0,30	pos	220,5
4	2 hour	0,20	pos	65,7
5	3 hour	0,17	pos	45,9
6	4 hour	0,18	pos	30,6
7	5 hour	0,20	pos	18,9
8	6 hour	0,24	pos	14,7
9	7 hour	0,25	pos	11,9
10	8 hour	0,28	pos	6,5

User 5 dated 25.06.2007				
Nr.	Time for collecting saliva after smoking	New Rapid STAT Collector Mirco flocced swab		
		quantity of saliva ml	result Rapid STAT THC	Result GC/MS THC ng/ml
2	30 min.	0,26	pos	780,5
3	1 hour	0,28	pos	345,8
4	2 hour	0,20	pos	110,9
5	3 hour	0,18	pos	39,6
6	4 hour	0,24	pos	50,6
7	5 hour	0,25	pos	31,3
8	6 hour	0,24	pos	34,9
9	7 hour	0,21	pos	20,5
10	8 hour	0,24	pos	9,8
11	9 hour	0,28	pos	8,5
12	10 hour	0,30	pos	4,7